ScaleReady.

G-Rex®: An Integrated Development and Manufacturing Platform for Cell and Gene-Modified Cell Therapy Applications

"Simplicity is the ultimate sophistication" - Leonardo Da Vinci

© 2020 Copyright ScaleReady

History of Cell Culture

An inverse relationship between oxygen availability and media volume (nutrient availability) necessitated constant interventions







Gas permeable bags corrected the inverse relationship, but continued to necessitate constant interventions





From VueLife Website

WHY ONE CENTIMETER: The selection of one centimeter thickness was based on research which indicated that human lymphocytes would expand at a maximal rate when the media fill was one centimeter thick. Greater thickness provided no greater expansion and is accordingly a waste of media. Lesser fill thickness resulted in lower expansion rates and is accordingly a waste of time.

Complex bioreactors are not a solution to solving scalable manufacturing because...













Development

- X No native small-scale model
- ! Unpredictable translation of small-scale data
- Excessive time and money spent attempting to create poorly predictive scale down models
- Excessive time and money spent optimizing CMC at full-scale

WHY REINVENT THE WHEEL WHEN YOU DON'T HAVE TO?



Manufacturing

- X Continuous fluidic connection of drug product
- ! Constraint is the manufacturing process itself
- ! Excessive investment in capital equipment for minimal throughput
- Limited potential for further innovation









- ! Overly complex movement of fluid
- Inefficient use of space
- No native small-scale model
- ! No high throughput manufacturing







- ! Overly complex movement of fluid
- ! Inefficient use of space
- ! No native small-scale model
- ! No high throughput manufacturing







- ! Overly complex movement of fluid
- ! Inefficient use of space
- No native small-scale model
- ! No high throughput manufacturing

Automated Instruments Sacrifice Scalability and Practicality

These needlessly complex instruments sacrifice scalability for automation and lack practicality when scaling out for larger patient populations.

Scale up

Small scale experiments don't easily scale up from benchtop to an instrument or "patient" dose



Scale out

Requiring one instrument per patient means large investment in cleanroom space & capital equipment



Strip away the needless complexity of a complex bioreactor and what do you get? A New Path Forward



G-Rex: The New Path Forward



Development

- ✓ Native small-scale model
- ! Predictable translation
- ! Linear path of least resistance to IND
- ! Cost- and time-efficient CMC improvements

Manufacturing

- ✓ De-couple drug product from instruments
- ! Constraints are individual unit operations
- ! Modest investment for high throughput
- ! Future-proofed



Surface Area





G-Rex[®] Theory

G-Rex[®]: How it Works

$G-Rex^{\otimes} = \underline{G}as$ Permeable <u>R</u>apid <u>Ex</u>pansion

Convection 10mL 1cm² Gas Permeable Membrane

Zone of Convection

Convection provides cells with nutrients on demand No mixing or perfusion needed Easily sample media for indicators of cell quantity

Zone of Diffusion

300 µm Boundary Layer of Diffusion Ideal for cell-to-cell communication

Allows cells to reside in a large volume of media to eliminate additional feeding

Unlimited Oxygen on Demand

G-Rex[®]: How it Works



G-Rex[®] has greatest range of cells/cm² (or cells/mL) in the industry



G-Rex[®]: How it Works



G-Rex® Best Practices

- Normalize to a Density (cells/cm²)
 - ✓ What governs production in G-Rex is space, not cell concentration.

Keep the cells on the G-Rex membrane

 Avoid resuspending the cells off the membrane, a microenvironment rich in oxygen & nutrients where the cells can secrete and receive the molecules they need to enter and sustain log-phase growth

> Use the full media volume capacity (for expansion)

- ✓ 10mL/cm² is the proper size nutrient pool and waste sink to facilitate maximum fold expansion without exchanging the media
- Trend Lactate
 - Lactate concentrations can be correlated to cell density, trended overtime, and used as a surrogate marker for cell proliferation

G-Rex[®]: An integrated development and manufacturing platform

G-Rex[®] Integrated Development and Manufacturing Platform for CGT

Linear Path of Least Resistance

- Continuity between R&D, PD, and GMP
- Optimize/troubleshoot quickly and cost efficiently with predictability

Modular Scale Out

- High-throughput manufacturing
 - ✓ 500,000 CAR-T doses in 70,000 sq.ft. facility
- Reduced investment in space, capital equipment, specialized labor, etc.





G-Rex "M" Series Integrated Development and Manufacturing Platform



ScaleReady Modular Manufacturing Platform



✓ Simplified high throughput manufacturing with a modest investment (low risk, high reward)

ScaleReady.

Proprietary & Confidential

Purpose Built Incubators for G-Rex500M-CS®







- Heavy duty shelving for standard incubators available today
- > $4x10^{11}$ cells in $<1m^2$

Description	120 V 50/60Hz	230V 50/60 Hz	100V (Japan only)
Heracell VIOS 250i CO ₂ Incubator			
Electropolished stainless steel interior with TC CO ₂ sensor	51033587	51033589	51033585
Electropolished stainless steel interior with IR CO ₂ sensor	51033597	51033599	51033605
Heracell Vios 250i CR (Cleanroom) CO ₂ Incubator CTS S	eries		
Electropolished stainless steel interior with IR CO ₂ sensor	51033783	51033784	51033782
Shelving System optimized for G-Rex [®] 500M-CS*			
VIOS 250i SST Shelving System to support 10 G-Rex [®] 500M		50164781	





Pilot Facility in Houston, TX

- Approximately 70,000 sq.ft.
- Standardized

Modular

Projected Annual Patient Capacity

- ➤ 30,000 TIL @ 80E9 cells/dose
- ➤ 500,000 CAR-T @ 4E9 cells/dose
- > 1,250,000 NK @ 2E9 cells/dose

G-Rex[®]: Optimized geometry enables simplified manufacturing

Simplifying Cell Culture with G-Rex[®]



✓ 10mL/cm² is the proper geometry to carry a culture from minimum to maximum density without the need to exchange the media



Frontload all 10mL/cm² to avoid resuspending cells and to maximize culture kinetics

Simplifying Cell Culture with G-Rex[®]





Vi = 100mL Protocol

Day 0: Stimulate 2E5 T cells and 2E7 feeder cells in 100mL

cytokine-supplemented media

Days 7, 8, 9, 10: Remove 100uL of media from top of well for lactate testing

Day 10: Re-suspend cells, remove 30uL cell suspension for counts

Vi = 16mL Protocol

Day 0: Stimulate 2E5 T cells and 2E7 feeder cells in 16mL cytokine-supplemented media

Day 4: Add 16mL cytokine-supplemented media

Day 7: Add 24mL cytokine-supplemented media, remove 30uL cell suspension for counts

Day 10: Add 44mL cytokine-supplemented media, remove 30uL cell suspension for counts



Stability of cytokines in media at different temperatures



GREX biotechne RD SYSTEMS

Experimental Workflow

Day 0

Thaw cells and

activate in G-Rex

*Sample media

Isolated T cells were thawed, plated, and activated in a G-Rex 6M at 0.5e6/mL/cm² using GMP Human T Cell Media supplemented with 5% hAB serum and 10 ng/mL IL-7 and 10 ng/mL IL-15. On Day 2, G-Rex was filled with fresh media to bring the total volume to 10 mL per cm² of G-Rex surface area. Media

Day 2

Fill G-Rex

with Media

*Sample media before

and after fill

 \rightarrow

Day 3-9

Fill and Forget Expand culture

without touching

*Sample media Day 5

Spike (Day 5) Expand culture with cytokine spike at Day 5 *Sample before and after spike

supernatant was collected for measuring cytokine concentrations on days 0, 2, and 9. Cytokine concentrations were quantified using the Ella Simple Plex automated immunoassay platform. All values shown were averaged across three donors.

Backend Analy

Ella quantification

média samples

Day 9

Cell counts and

phenotyping

*Sample media

Materials and Methods

Product	Material	Part #
G-Rex 6M plate*	Culture Vessel	80660M (ScaleReady)
GMP Human T Cell Media	Media	CCM038-GMP
GMP IL-7 (10 ng)**	Cytokine	BT-007-GMP
GMP IL-15 (10 ng)**	Cytokine	BT-015-GMP
Ella	Analytical Instrument	600-100
IL-7 (72x1)	Simple Plex Cartridge	SPCKB-PS-000506
IL-15 (72x1)	Simple Plex Cartridge	SPCKB-PS-000500

Available exclusively from ScaleReady

** Available from ScaleReady G-Rex is a registered trademark of Wilson Wolf Manufacturing, LLC





 G-Rex[®] process can be simplified using R&D Systems reagents to a single bolus on D0 to achieve max cell density with desirable phenotypes



Fill & Forget with G-Rex[®]

Fill-and-Forget Workflow Produced Equivalent T Cell Expansion and Viability



Both workflows supported cell expansion to reach confluent

No cytokine control showed minimal expansion as expected.

density of 30-40 million cell/cm².

Data labels denote average viability.

Average Cell Density at Harvest

Fill-and-Forget Workflow Produced Equivalent Expanded T Cell Phenotypes



IL-7/15 Average Phenotype (CCR7 vs CD45RA)

G-Rex[®]: Linear scalability enables continuity

Linear Scalability of G-Rex[®] (k562)







- Linear scalability is unique to the G-Rex platform.
- ✓ Streamline and shorten the lifecycle of a cell therapy by seamlessly scaling up from RUO well plates to a GMP bioreactor
- ✓ Take an established process back to the drawing board by scaling down from your GMP production vessel to the G-Rex6M well plate.





Proprietary & Confidential 38

*Unpublished data from Dr. Eric Tran Lab. Shared under Data Use Agreement



CD4+ T cells







Linear Scalability of G-Rex[®] (TCR-T)



Linear scalability of G-Rex[®] enables process continuity


Linear scalability of G-Rex[®] enables product continuity



Linear scalability of G-Rex[®] enables product continuity





Process Continuity = **Product** Continuity



Figure 1. NK-cell expansion rate is comparable between Grex 6M, 100M, and 500M vessels.

Figure 2. NK-cell potency is comparable between Grex 6M, 100M, and 500M vessels.

Linear Scalability of G-Rex[®] (NK)





✓ Comparable Expansion Kinetics, Post Thaw Viability and Anti-Tumor Activity of NK Cells When Cultured in a G-Rex6M or G-Rex500M

ScaleReady.

Proprietary & Confidential 45

Linear Scalability of G-Rex[®] (T cell)



✓ Comparable Expansion Kinetics, Viability and Antigen Specificity of T Cells When Cultured in a G-Rex10M, G-Rex100M or G-Rex500M

ScaleReady,

Proprietary & Confidential 46

Cell)Readv







✓ Despite donor-to-donor variability, the G-Rex[®] is remarkably scalable from 6M up to 500M



Days expansion

Linear Scalability of G-Rex[®] (CAR-T)



G-Rex6M-2cm² vs G-Rex6M



 Across 3 different healthy donors, lentivirus transduced T cells display similar growth kinetics in the G-Rex6M-2cm² and G-Rex6M. Similar cell densities (cells/cm²) were reached upon harvest.

The G-Rex6M-2cm² shows comparability to the G-Rex6M – but requires only 20% of the reagents

Linear Scalability of G-Rex[®] (CAR-T)





- ✓ Upscaling the process from the G-Rex6M to manufacturing scale in the G-Rex100M does not impact transduction efficiencies, cell yields, or purity
- Process development in the G-Rex6M can be translated directly to the G-Rex100M for product manufacture

() **周** TakaRa **Cell**Ready



✓ FMC63 anti-CD19 CAR-T manufacturing scalability using retrovirus and RetroNectin-coated G-Rex[®]

✓ Cell Ready compared transduction/expansion in the G-Rex6M to the G-Rex100M (n = 1 donor)

✓ 2-day activated PBMCs transduced overnight (0.5x10⁶ cells/cm²), expanded until Day 8 for harvest

G-Rex[®]: Flexibility or and Consistency







Figure 2 Identifying the optimal seeding density to support maximum cell output. Panel (a) shows the expansion of K562 cells in G-Rex devices that were initiated with different seeding densities (0.0025, 0.125, 0.25, 0.50, and 1.0×10^6 cells/cm²). A half medium change was performed every day in all conditions. Panel (b) shows the final cell number on day 14 of culture (reported as cells/cm²). Panel (c) shows the fold increase in the cell numbers on day 9.

Flexible Seeding Densities (input)

Consistent Max Cell Densities (output)

Flexibility & Consistency with G-Rex[®] (CAR-T)



GREX Cityof Hope.

NOVEL, EFFICIENT SYSTEM FOR CULTURE OF TRANSDUCED HEMATOPOIETIC STEM PROGENITOR CELLS (HSPC) FOR VECTOR COPY NUMBER DETERMINATION

Monica Torres-Coronado¹, Ilse Arciga¹, Diana L. Browning¹, Xiu-Li Li¹, Agnes Gardner¹, Sergio Branciamore², John A. Zaia¹, Angelo A. Cardoso¹

Center for Gene Therapy ¹ and Department of Computational and Quantitative Medicine ² of City of Hope, Duarte, CA

Establishing Seeding Conditions to culture HSPC in G-Rex Plates (TDX with GFP LVV)





55



G-Rex[®]: Trending Metabolites

Measuring Lactate with G-Rex[®]





- ➤ Lactate concentrations in G-Rex[®]...
 - are useful surrogate measurements and leading indicators
 - correlates nearly perfectly to cell density
 - can inform harvest or passage events

ScaleReady.

Gagliardi C, Khalil M, Foster AE. Streamlined production of genetically modified T cells with activation, transduction and expansion in closed-system G-Rex bioreactors. Cytotherapy. 2019 Dec;21(12):1246-1257. DOI: 10.1016/j.jcyt.2019.10.006.

G-Rex[®]: Closed System Overview







- 1 Gas Permeable Membrane
- 2 Polycarbonate Shell (Vessel)
- 3 Polypropylene Multi-Port Cap
- 4 Vent Filter (0.2 µm)
- 5 Sample Line Tubing
- 6 MicroClave® Connector
- 7 Reduction Line Tubing
- 8 PVC Weldable Reduction Line
- 9 Harvest Line Tubing
- 10 PVC Weldable Harvest Line

GatheRex[®] for efficient, closed system G-Rex[®] processing





Semi-automated volume reduction for concentrated cell product

Semi-automated cell harvest with high recovery



Figure 7 Automated collection of cells in a closed G-Rex system. Panel (**a**) shows the setup of "GatheRex" draining the excess medium to concentrate the cells: (1) pump, (2) 0.2-µm sterile filter, (3) medium collection button, (4) medium line clamp, (5,6) medium/supplements input ports, (7) medium sampling port, (8) medium collection port, (9) medium bag, (10) medium, and (11) medium line optical detector. Panel (**b**) shows the GatheRex harvesting the concentrated cells: (12) resuspended cells, (13) cell collection button, (14) cell line clamp, (15) cell collection port, (16) cell collection bag, (17) harvested cells, and (18) cell line optical detector. Panel (**c**) shows a comparison in cell recovery achieved when cells are collected manually versus by the GatheRex. To confirm the absence of cells in cell collection tubing, the medium collection bag or the G-Rex following cell collection with the GatheRex, each component was washed with media and residual cells were harvested and counted.



Grex500 Volume Reduction





G-Rex500M-CS



- ✓ Use large diameter C-Flex tubing to quickly volume reduce a G-Rex500M-CS or G-Rex500M-TF
 - Requires GatheRex model 80000Z
- Volume reduce 90% (4.5L) of waste media out in less than 3 minutes
 - ✓ ~80% faster when compared to small diameter PVC tubing
- Concentrate cells to 40e6 cells/mL in ~500mL for faster downstream processing

G-Rex[®] for CAR-T cell therapy



G-Rex[®] for Autologous CAR-T

ScaleReady.





Fill up the G-Rex with the full volume of media (10mL/cm²) and go from *minimum cell density* to *maximum cell density* in a single REP with minimal (potentially no) intervention

A single G-Rex100M-CS (below) will produce
3+ billion cells

A standard, stackable incubator (like Thermo Fisher's Heracell250i) can hold ~45 G-Rex100M-CS devices per incubator



Bajgain P, Mucharla R, Wilson J, et al. Optimizing the production of suspension cells using the G-Rex "M" series. *Mol Ther Methods Clin Dev.* 2014;1:14015. Published 2014 May 14. doi:10.1038/mtm.2014.15





G-REX[®] centric, cGMP Compliant, Non-Viral CAR-T cell manufacturing

ScaleReady,

Shy, B. R., Vykunta, V. S., Ha, A., Talbot, A., Roth, T. L., Nguyen, D. N., Pfeifer, W. G., Chen, Y. Y., Blaeschke, F., Shifrut, E., Vedova, S., Mamedov, M. R., Chung, J.-Y. J., Li, H., Yu, R., Wu, D., Wolf, J., Martin, T. G., Castro, C. E., ... Marson, A. (2022). High-yield genome engineering in primary cells using a hybrid ssdna repair template and small-molecule cocktails. *Nature Biotechnology*, *41*(4), 521–531. https://doi.org/10.1038/s41587-022-01418-8

G-Rex[®] for Autologous CAR-T





G-REX[®] delivers enough CAR⁺ cells for patient dose

G-REX[®] delivers high quality cells with favorable phenotypes

Shy, B. R., Vykunta, V. S., Ha, A., Talbot, A., Roth, T. L., Nguyen, D. N., Pfeifer, W. G., Chen, Y. Y., Blaeschke, F., Shifrut, E., Vedova, S., Mamedov, M. R., Chung, J.-Y. J., Li, H., Yu, R., Wu, D., Wolf, J., Martin, T. G., Castro, C. E., ... Marson, A. (2022). High-yield genome engineering in primary cells using a hybrid ssdna repair template and small-molecule cocktails. *Nature Biotechnology*, *41*(4), 521–531. https://doi.org/10.1038/s41587-022-01418-8

ScaleReady.

Proprietary & Confidential 66





- G-REX[®] is utilized extensively in academic research institutions to streamline the (linear) pathway from benchtop to bedside
- G-REX[®] is ideal for viral and nonviral CAR-T therapies

ScaleReady.

Balke-Want, H., Keerthi, V., Gkitsas, N., Mancini, A., Kurgan, G., Fowler, C., Xu, P., Liu, X., Asano, K., Patel, S., Fisher, C., Brown, A., Tunuguntla, R., Patel, S., Sotillo, E., Mackall, C., & Feldman, S. (2023). Gene editing/gene therapies: Homology-independent targeted insertion (HITI) enables guided car knock-in and efficient clinical scale CAR-T cell manufacturing. Cytotherapy, 25(6). <u>https://doi.org/10.1016/s1465-3249(23)00166-4</u>

Proprietary & Confidential 67

G-Rex[®] for Autologous CAR-T





centrifuge. The use of a large format G-Rex device allows for all major steps to be carried out on the day of initiation (Figure 1 and Figure 6), and only sterility and product qualification steps near the end of the culture period would require further technical handling. While a clinical protocol incorporating the advances presented here is yet to implemented in a clinical trial, we present here all the steps required to build a low-cost CAR-T pointof-care production pipeline, at scale, that can be implemented in a global setting either in a dedicated clean-room area or in a mobile cleanroom setting.



Cell Viability











ScaleReady,

Xiong, Y., Xie, Y., Zhu, Z., Oparaocha, I., Sleesareva, O., Dropulić, B., & Orentas, R. (2023). Technological and Manufacturing Innovation Drive improved access to engineered T cell therapies. Medical Research Archives, 11(11). <u>https://doi.org/10.18103/mra.v11i11.4759</u>

G-Rex[®] for Autologous CAR-T



- CAR-T manufactured in G-REX[®] results in high expansion of functional CAR-T cells.
- "The use of G-REX[®] resulted in high expansion of CAR-T cells (>180 fold) with a high number of central memory phenotypes that were shown to be critical for CAR-T cell persistence and increased survival."



ScaleReady.

T cells

HeLa

120,000-

90,000

60,000

30,000

Our G-Rex[®] CAR-T platform that is ready-to-scale



- TcBuster reagents are introduced into the cells via 1. electroporation
 - Transposase is mRNA to limit DNA toxicity and ensure safety •
 - Transposon is the Nanoplasmid[™] vector with minimal • backbone to improve transposition rates and eliminate encoded antibiotic marker
- TcBuster mRNA is translated into transposase enzyme 2.
- TcBuster transposase binds to ITR regions on transposon 3. Nanoplasmid[™] vector
- TcBuster cuts gene of interest (GOI) out of transposon 4.
- TcBuster inserts GOI into genome 5.
- GOI transcribed and translated 6.
- GOI stably expressed in cells

ScaleReady.



✓ Adjusted Copy Number
✓ Potency

- Sacrificial wells in G-Rex 6M well plate
- Leaves the G-Rex 100M left untouched



Viability Monitoring





ScaleReady.

Proprietary & Confidential



Test (%)	Donor 1 Results	Donor 2 Results	Donor 3 Results
Post-thaw Viability	93.8	80.2	96.9
Cell Count (% recovery at thaw)	70	73	99



Test (%)	Donor 1 Results	Donor 2 Results	Donor 3 Results
T cell Purity (CD4 + CD8)	97.9%	99.0%	96.3%

ScaleReady.

Proprietary & Confidential



Bolded Values used for T cell potency assay normalization and copy number assay

ScaleReady.

Proprietary & Confidential

Transposon copy number determination by dPCR



Test (%)	Donor 1 Results	Donor 2 Results	Donor 3 Results
Copy Number (Adjusted to CAR+)	5.6	5.5	5.8
Non-viral Manufacturing of CAR T cells in G-Rex®



Non-viral Manufacturing of CAR T cells in G-Rex®



Transposon Batch Size	Patients Treated at 50 ug/mL EP	Patients Treated at 150 ug/mL EP
100 mg	400	133
250 mg	1,000	333
500 mg	2,000	666
1 g	4,000	1,333

*Assumes 2E8 cells/mL and 1E9 total cells processed (5 x 1 mL electroporation reactions)

Non-viral Manufacturing of CAR T cells in G-Rex®

TcBuster-M[™] Advantages

Larger cargo delivery allows for more efficacious products
 Improved timelines for discovery and manufacturing
 Integrates genes in safer locations than lentiviral methods
 Scalable from research to commercial scale
 Reduced COGS versus viral editing
 Stronger supply chain for manufacturing of nucleic acids

G-Rex[®] for Tumor Infiltrating Lymphocyte (TIL) therapy

G-Rex[®] in world's first TIL Therapy





Dr. Steven Rosenberg in his National Cancer Institute lab with the legacy G-Rex100

Simplified Method of the Growth of Human Tumor Infiltrating Lymphocytes (TIL) in Gas-Permeable Flasks to Numbers Needed for Patient Treatment

Jianjian Jin¹, Marianna Sabatino¹, Robert Somerville², John R. Wilson³, Mark E. Dudley², David F. Stroncek¹, and Steven A. Rosenberg²

¹Cell Processing Section, Department of Transfusion Medicine, Clinical Center, National Institutes of Health, Bethesda, Maryland, USA

²Surgery Branch, National Cancer Institute, National Institutes of Health, Bethesda, Maryland, USA Immunotherapy

³Wilson Wolf Manufacturing, New Brighton, MN, USA



- ✓ Far fewer vessels are needed because G-Rex can support much greater cell densities/concentrations
- \checkmark 3-4 fold reduction in the use of media, cytokines, activating reagents, serum, etc.
- Reduced culture/harvest volume reduces need for specialized cell process equipment
- Less incubator space required \checkmark
- \checkmark Less labor is required because laboratory staff maintain fewer vessels, prepare less media, vessels do not require media exchanges, less volume is processed post REP

ScaleReady.

Jin, J., Sabatino, M., Somerville, R., Wilson, J. R., Dudley, M. E., Stroncek, D. F., & Rosenberg, S. A. (2012). Simplified method of the growth of human tumor infiltrating lymphocytes in gas-permeable flasks to numbers needed for patient treatment. Journal of Immunotherapy, 35(3), 283-292. https://doi.org/10.1097/cji.0b013e31824e801f





G-Rex100M-CS



G-Rex500M-CS





Figure 1

(57) Abstract: The present invention provides improved and/or shortened methods for expanding TILs and producing therapeutic populations of TILs, including novel methods for expanding TIL populations in a closed system that lead to improved efficacy, improved phenotype, and increased metabolic health of the TILs in a shorter time period, while allowing for reduced microbial contamination as well as decreased costs. Such TILs find use in therapeutic treatment regimens.

✓ G-Rex[®] is a core technology in the world's first approved TIL therapy

✓ G-Rex[®] is the gold standard for TIL therapy



The beneficial effects of a gas-permeable flask for expansion of Tumor-Infiltrating lymphocytes as reflected in their mitochondrial function and respiration capacity

Marie-Andrée Forget^{1,†}, Cara Haymaker^{1,†}, Jennifer B Dennison², Christopher Toth¹, Sourindra Maiti³, Orenthial J Fulbright¹, Laurence J N Cooper³, Patrick Hwu¹, Laszlo G Radvanyi^{1,4,5}, and Chantale Bernatchee^{1,*}

¹Department of Melanoma Medical Oncology; The University of Texas MD Anderson Cancer Center (MDACC); Houston, TX USA; ²Department of Systems Biology; The University of Texas MDACC; Houston, TX USA; ⁴Lion Biotechnologies; Tampa, FL USA; ⁵Department of Immunology; H. Lee Moffitt Cancer Center; Tampa, FL USA

Interestingly, we observed significantly better growth when using the G-Rex flask for TIL lines that did poorly (lower range) in the traditional devices (Fig. 1B, n = 4, p = 0.007). This observation explains the "trend" of better growth observed in the G-Rex flask across all samples since good growing TIL lines were not affected by the type of device used (Fig. 1). Neither device negatively affected the post-REP viability (>85 %, data not shown). The differential growth of poor growing TIL lines in the two systems used prompted us to explore potential underlying explanations.



Figure S1. TIL fold expansion on day 14 of the Rapid Expansion Protocol. TIL fold expansion obtained for the first 45 patients treated with TIL ACT at MDACC. TIL were expanded in T175 flasks for the first week and transferred to 3L cell culture bags for the last 7 days. The mean fold expansion was 1651 with a standard deviation of 799. Poor growth was established as TIL expansion not reaching the mean minus one standard deviation (1651 - 799 = 852).

ScaleReady.

Forget, M.-A., Haymaker, C., Dennison, J. B., Toth, C., Maiti, S., Fulbright, O. J., Cooper, L. J., Hwu, P., Radvanyi, L. G., & Bernatchez, C. (2015). The beneficial effects of a gas-permeable flask for expansion of tumor-infiltrating lymphocytes as reflected in their mitochondrial function and respiration capacity. *Oncolmmunology*, *5*(2). https://doi.org/10.1080/2162402x.2015.1057386

G-Rex[®] Preserves Clonal Diversity





Figure 2. Rapid expansion of TIL in the Gas-permeable flask (G-Rex) does not favor selection and expansion of specific T cell clones. Clonal diversity was evaluated by measurement of the expression of major TCR α and β chain gene by the NanoString nCounter[®] technology. The analysis of the level of expression of 45 V α and 46 V β genes in RNA isolated from the pre-REP TIL lines in comparison with the different expansion devices demonstrated no difference in the clonal diversity obtained post REP. Two out of 4 TIL lines from 4 melanoma patients are shown. Statistical analysis using a Spearman correlation comparing cells grown in traditional flask and bag (red) vs. Gas-permeable flask (blue) is shown for both V α and V β genes.

ScaleReady,

Forget, M.-A., Haymaker, C., Dennison, J. B., Toth, C., Maiti, S., Fulbright, O. J., Cooper, L. J., Hwu, P., Radvanyi, L. G., & Bernatchez, C. (2015). The beneficial effects of a gas-permeable flask for expansion of tumor-infiltrating lymphocytes as reflected in their mitochondrial function and respiration capacity. *Oncolmmunology*, *5*(2). https://doi.org/10.1080/2162402x.2015.1057386

The utility of G-Rex[®] on TIL cell quality





Figure 3. Bioenergetic analysis of melanoma TIL lines to evaluate mitochondrial function. (**A–D**) OCR of 4 TIL lines on day 7 of the REP (**A and B**) or sorted, pre-REP CD8⁺BTLA⁺ and CD8⁺BTLA⁻ TIL line from 3 patients (**C and D**) were determined using a Seahorse XP96 Bioanalyzer. OCR were calculated after 3 min of mix time and 4 min of measurement time. OCR[Mito] is the total mitochondrial OCR, the value just prior to the Oligomycin (OLG) injection minus the non-mitochonddrial OCR component determine by Antimycin (AA) treatment. OCR[OLG] is the component of the OCR that is sensitive to Oligomycin treatment, the rate used by ATP-syntase. The OCR on a per cell basis was determined by dividing the OCR by the seeded TIL count of 250000. (**E**) OCR/ECAR ratio is the OCR[OLG] divided by the ECAR basal which is reflextive of a cell dependence on glycolysis. (**F**) pre-REP TIL were stained using a MitoTracker dye. The MFI of the MitoTracker staining of the CD8⁺BTLA⁺ and CD8⁺BTLA⁻ subset from 3 TIL lines is shown. The histogram shows staining of a representative TIL line. Statistical significance was determined by using a paired *t*-test.

- ✓ The G-Rex[®] grown cells consistently had a higher respiratory capacity than the corresponding T175 flask cells. All together, the metabolic phenotype of melanoma TIL propagated in G-Rex[®] was consistently more oxidative.
- High spare respiratory capacity is believed to be a feature of memory T cells that imparts them with **better in vivo persistence**, as opposed to effector T cells. We are thus postulating that the high spare respiratory capacity and the high OCR that we observed in TIL grown in high oxygen environment, more specifically the CD8+BTLA+ cells, could therefore result in **enhanced survival and persistence** after transfer, which could confer **better tumor control.**

ScaleReady.

Forget, M.-A., Haymaker, C., Dennison, J. B., Toth, C., Maiti, S., Fulbright, O. J., Cooper, L. J., Hwu, P., Radvanyi, L. G., & Bernatchez, C. (2015). The beneficial effects of a gas-permeable flask for expansion of tumor-infiltrating lymphocytes as reflected in their mitochondrial function and respiration capacity. *Oncolmmunology*, *5*(2). https://doi.org/10.1080/2162402x.2015.1057386

G-Rex[®] Simplifies TIL Manufacturing





unavailable. In summary, in addition to production of greater numbers of viable TILs, transition to G-Rex **100MCS** flasks has resulted in fewer total vessels, less media, less incubator space and less labor than traditional culture methods using flasks and bags.

Given the benefits shown by addition of anti-4-1BB antibody during preREP and culturing TILs in G-Rex 100MCS during REP, we designed a phase 1 clinical trial to incorporate both process changes. When taken together, the optimization decreased the culture period from 51-37 days and limited the technical intervention to 14 days. Most importantly, optimizing the TIL culture expansion allowed for consistent achievement of the optimal therapeutic dose of 60×10^9 cells for all patients compared with achievement of optimal dose in only 40% of patients using traditional culture methods. This clinical trial









Proprietary & Confidential 87

Hopewell, E. L., Cox, C., Pilon-Thomas, S., & Kelley, L. L. (2019). Tumor-infiltrating lymphocytes: Streamlining a complex manufacturing process. Cytotherapy, 21(3), 307-314. https://doi.org/10.1016/j.jcyt.2018.11.004

G-Rex[®] for TCR-T cell therapy

ScaleReady,

G-Rex[®] for Autologous TCR-T







Central memory T cells seem to be required for longterm in vivo repopulation, and there are concerns that excessive in vitro T-cell proliferation before infusion may lead to terminal differentiation and exhaustion. To determine whether the G-Rex would favor the unwanted production of exhausted, terminally differentiated effector cells,⁴⁰ we measured both cell proliferation and cell death. We found that the increased cell numbers were a result of improved survival rather than increased proliferation. In confirmation, phenotypic analysis of G-Rex-grown CTL revealed no apparent detrimental effects on the ratio of effector:memory phenotype.

- ✓ A single G-Rex500M-CS (above) produces 15-20 billion cells
- ✓ Standard, stackable incubators hold ~8 devices. This equals 240-320 billion cells in only two stacked incubators!
- High viability rates means less population doublings are required to reach final cell numbers
- Less population doublings yield move favorable, less differentiated, phenotypes

Vera JF, Brenner LJ, Gerdemann U, Ngo MC, Sili U, Liu H, Wilson J, Dotti G, Heslop HE, Leen AM, Rooney CM. Accelerated production of antigen-specific T cells for preclinical and clinical applications using gas-permeable rapid expansion cultureware (G-Rex). J Immunother. 2010 Apr;33(3):305-15. doi: 10.1097/CJI.0b013e3181c0c3cb. PMID: 20445351; PMCID: PMC2946348.

Jin, J., Gkitsas, N., Fellowes, V.S. *et al.* Enhanced clinical-scale manufacturing of TCR transduced T-cells using closed culture system modules. *J Transl Med* **16**, 13 (2018). https://doi.org/10.1186/s12967-018-1384-z



90





Linear Scalability of G-Rex[®] (TCR-T)



Simplifying Cell Culture with G-Rex[®]





Vi = 100mL Protocol

Day 0: Stimulate 2E5 T cells and 2E7 feeder cells in 100mL

cytokine-supplemented media

Days 7, 8, 9, 10: Remove 100uL of media from top of well for lactate testing

Day 10: Re-suspend cells, remove 30uL cell suspension for counts

Vi = 16mL Protocol

Day 0: Stimulate 2E5 T cells and 2E7 feeder cells in 16mL cytokine-supplemented media

Day 4: Add 16mL cytokine-supplemented media

Day 7: Add 24mL cytokine-supplemented media, remove 30uL cell suspension for counts

Day 10: Add 44mL cytokine-supplemented media, remove 30uL cell suspension for counts

ScaleReady.

G-Rex[®] for γδ T cell therapy



Luminary Therapeutics' Gamma 2.0 Manufacturing achieves Enhanced Expansion





Client Testimony

- Luminary is utilizing the G-Rex500M-CS and starting from ~20% of a full leukopak has achieved 194 patient doses of 65% CAR+ engineered gamma delta cells
- Patient Dose size calculated at 350M cells
- Luminary is driving new levels of patient access in cell therapy



- "Our manufacturing platform when coupled with the G-Rex will allow us to make
 50,000 patient products per year in two conventional clean rooms" Luminary CEO
- "We are confident that we are currently operating with the lowest COGS for any cell therapy product in the entire industry" – Luminary CEO

G-Rex[®] for Natural Killer (NK) cell therapy

ScaleReady,

Proprietary & Confidential 96

G-Rex[®] for Donor-Derived Allo NK





Cytotherapy. 2012 October; 14(9): 1131-1143. doi:10.3109/14653249.2012.700767.

Large-scale *ex vivo* expansion and characterization of natural killer cells for clinical applications

NATALIA LAPTEVA¹, APRIL G. DURETT¹, JIALI SUN¹, LISA A. ROLLINS¹, LESLIE L. HUYE¹, JIAN FANG¹, VARADA DANDEKAR¹, ZHUYONG MEI¹, KIMBERLEY JACKSON¹, JUAN VERA¹, JUN ANDO¹, MINHTRAN C. NGO¹, ELAINE COUSTAN-SMITH³, DARIO CAMPANA³, SUSANN SZMANIA⁴, TARUN GARG⁴, AMBERLY MORENO-BOST⁴, FRITS VANRHEE⁴, ADRIAN P. GEE^{1,2}, and CLIONA M. ROONEY^{1,2}

 $^1 \text{Center}$ for Cell and Gene Therapy, The Methodist Hospital, Texas Children's Hospital, Houston, TX

²Department of Pediatrics, Baylor College of Medicine, Houston, Texas, USA

³National University of Singapore, Singapore

⁴University of Arkansas Medical Center, Little Rock, Arkansas, USA

Figure 2.

G-Rex yields higher numbers of NK cells and requires fewer culture manipulations than

eultures in bags. (A) Schematic view of NK cell expansion in bags and G-Rexes. Less cell processing is required to grow a similar number of cells in G-Rexes. (B) NK cells from three different donors were seeded in the 197-mL bags and G-Rex100s. Bags were seeded with PBMC containing 2×10^6 NK cells and 2×10^7 K562-mb15-41BBL cells in a 40-mL starting volume and cultures required feedings with 40 mL fresh media and cytokine every other day, i.e. days 2, 4, 6 and 8, resulting in 200 mL of total volume by day 8. G-Rex100s were seeded with PBMC containing 2×10^6 NK cells and 2×10^7 K562-mb1L15-41BBL in 400 mL medium and were not subsequently fed. Cells were counted and analyzed by flow cytometry for percentage CD56⁺ CD3⁻ cells (n = 3). (C) Fold expansion of CD56⁺ CD3⁻ cells in G-Rexes and bags (n = 3). *P < 0.5.

G-Rex[®] for Umbilical Cord Blood derived NK



B 10000 STORE ALLS Fresh + IL-2 Fresh + aAPC Frozen + aAPC Frozen + aAPC 1000 100 100 0.1 0.1 0.01 Days



Figure 2. Co-culture of CB MNCs with IL-2 and aAPCs yields significantly greater expansion of NK cells than culture with IL-2 alone. A. Mean fold growth of CD56⁺/CD3⁻ NK cells from 8 fresh and 6 frozen cord blood expansions with aAPCs and IL-2 versus 3 expansions with IL-2 alone (14 day culture). B. Time course of NK cell growth over 14 day culture between all 3 conditions. By day 7, the fresh CB aAPC-containing culture demonstrated greater NK cell growth than culture with IL-2 alone (p<0.05). The frozen CB showed a similar trend at day 7, which did not reach statistical significance (p=0.06). C. All three culture conditions yielded comparable, low percentages of CD3⁺ cells: 0.44%, 0.74% and 0.66% CD3⁺ cells from the culture with IL-2 alone, fresh CB MNCs with aAPC feeders or frozen CB MNCs with aAPC feeders respectively (p>0.5 for all comparisons). Mean +/- SD is shown for each figure. P<0.05 where indicated (*). doi:10.1371/journal.pone.0076781.g002

Antigen Presenting Cell-Mediated Expansion of Human Umbilical Cord Blood Yields Log-Scale Expansion of Natural Killer Cells with Anti-Myeloma Activity

Nina Shah¹*, Beatriz Martin-Antonio¹, Hong Yang¹, Stephanie Ku², Dean A. Lee³, Laurence J. N. Cooper³, William K. Decker^{2,4}, Sufang Li¹, Simon N. Robinson¹, Takuya Sekine¹, Simrit Parmar¹, John Gribben⁵, Michael Wang⁶, Katy Rezvani¹, Eric Yvon¹, Amer Najjar⁷, Jared Burks⁸, Indreshpal Kaur¹, Richard E. Champlin¹, Catherine M. Bollard², Elizabeth J. Shpall¹

1 Department of Stem Cell Transplantation and Cellular Therapy, The University of Texas M.D. Anderson Cancer Center, Houston, Texas, United States of America, 2 Center for Cell and Gene Therapy, Baylor College of Medicine, Houston, Texas, United States of America, 3 Department of Pediatrics, The University of Texas M.D. Anderson Cancer Center, Houston, Texas, United States of America, 4 Department of Pathology and Immunology, Baylor College of Medicine, Houston, Texas, United States of America, 5 Institute of Cancer, Queen Mary University of London, Centre for Medical Oncology, Barts and The London School of Medicine, London, United Kingdom, 6 Department of Lymphoma, The University of Texas M.D. Anderson Cancer Center, Houston, Texas, United States of America, 7 Department of Experimental Diagnostic Imaging, The University of Texas M.D. Anderson Cancer Center, Houston, Texas, United States of America, 8 Department of Leukemia Research, The University of Texas M.D. Anderson Cancer Center, Houston, Texas, United States of America

ScaleReady,



Making Cancer History®

G-Rex[®] vs. WAVE for NK cells





FIG. 4: Comparison of NK cell expansion in G-Rex flasks and WAVE bioreactor. PBMCs from healthy donor containing 5×10⁶ NK cells were co-cultured with K562-41BBL-mbIL-15 cells for 9 days. Cells were moved to WAVE on day 4 of culture. (A, B, C) Similar numbers and fold of NK expansion were observed in G-Rex and WAVE cultures. (D) Potency (killing K562 cells) of both products was also comparable. Similar expansion was obtained for three donors and the representative data for one donor is shown. Critical ReviewsTM in Oncogenesis, 19(1-2):121–132 (2014)

Clinical Grade Purification and Expansion of Natural Killer Cells

Natalia Lapteva^{1*}, Susann M. Szmania², Frits van Rhee², & Cliona M. Rooney^{1*}

¹ Center for Cell and Gene Therapy, Baylor College of Medicine, Houston, TX 77030; ²University of Arkansas Medical Center, Little Rock, AK 72205

toxicity assays). Interestingly, the WAVE-generated products contain fewer CD3⁺ T cells and a higher frequency of CD56⁺CD3⁻ NK cells, perhaps because T cells prefer static culture (Fig. 4). The clinical trial in multiple myeloma demands both autologous and allogeneic NK cells. In the case of allogeneic products, CD3⁺ T cells are efficiently removed (as well as CD3+56+ cells) by depletion with the CliniMACS CD3 reagent, reducing the CD3+ cell frequency from 19±9% to 0.4±0.5% (N=5). This depletion allows less than 5×10⁵ CD3⁺ CD56⁻ T cells kg⁻¹ for infusion, as mandated by the clinical protocol.

ScaleReady.





Supplemental Figure 5. A. Fold expansion of NKF-NK cells at 5:1 feeder-to-NK ratio and 200U/mL IL-2 after 5 weeks, n=6. B, Expansion of NK cells using NKF cells at 5:1 feeder-to-NK ratio and 200U/mL IL-2 for 2 weeks, using G-Rex flasks, n=8. ***p<0.001. Data represent mean +/- SEM.

✓ G-Rex[®] enables robust fold expansion of Natural Killer cells

ScaleReady,

G-Rex[®] for Donor-Derived Allo NK





 G-Rex[®] can be used with multiple feeder cell lines to produce large amounts of NK cells with high tumor cytotoxicity

Supplemental Figure 2. A. Fold expansion of K562-NK and NKF-NK cells at 5:1 feeder-to-NK ratio and 100U/ mL IL-2 after 2 weeks, n=1. B, The cytotoxic activity of 2 weeks-expanded NKF-NK and K562-NK was assessed against OCI-AML3 and HCT116 cells after 4hr co-culture, n=3. The NK cell-to-target ratio was 1-1. Data represent mean +/- SEM.

ScaleReady.

А





ScaleReady.

GMP-Compliant Universal Antigen Presenting Cells (uAPC) Promote the Metabolic Fitness and Antitumor Activity of Armored Cord Blood CAR-NK Cells

Enli Liu^{1†}, Sonny O. T. Ang^{1†}, Lucila Kerbauy^{1,2,3}, Rafet Basar¹, Indreshpal Kaur¹, Mecit Kaplan¹, Li Li¹, Yijiu Tong¹, May Daher¹, Emily L. Ensley¹, Nadima Uprety¹, Ana Karen Nunez Cortes¹, Ryan Z. Yang¹, Ye Li¹, Hila Shaim¹, Francia Reyes Silva¹, Paul Lin¹, Vakul Mohanty⁴, Sunil Acharya¹, Mayra Shanley¹, Luis Muniz-Feliciano¹, Pinaki P. Banerjee¹, Ken Chen⁴, Richard E. Champlin¹, Elizabeth J. Shpall¹ and Katayoun Rezvani^{1*}

¹ Department of Stem Cell Transplantation and Cellular Therapy, The University of Texas MD Anderson Cancer Center, Houston, TX, United States, ² Departments of Stem Cell Transplantation and Hemotherapy/Cellular Therapy, Hospital Israelita Albert Einstein, Saő Paulo, Brazil, ³ Department of Genetics and Evolutionary Biology, Human Genome and Stem Cell Research Center, Biosciences Institute, University of Saő Paulo, Saő Paulo, Brazil, ⁴ Department of Bioinformatics and Computational Biology, The University of Texas MD Anderson Cancer Center, Houston, TX, United States

GMP NK Cell Expansion With uAPC

NK cells were expanded using the G-Rex[®] bioreactor (Saint Paul, Minnesota, USA) in co-cultures with uAPC for the CB-NK CAR clinical trial (21) (NCT Number: NCT03056339). This closed cell culture system enables using flasks with a gas-permeable membrane, allowing for optimal gaseous exchange and ensuring aerobic growth kinetics. Using previously optimized variables which included cell seeding density, calibrated media volume, and media formulation, we were able to adapt the G-Rex M100 series from pre-clinical protocols to translate and linearly scale the procedure to the desired yield in the GMP setting.



CB derived, optimized affinity CD38 CAR-NK cells with CD38 KO show promising *in-vivo* activity in a Multiple Myeloma model

S Brophy¹, L Kirkham McCarthy¹, M Köylijärvi¹, D Hardwicke¹, H Hintz², N Otto², B Ettestad², D Hermanson², E Shevlin¹, M Reilly¹ and M O'Dwyer¹ Bio-Techne² and ONK Therapeutics¹



- ✓ The process resulted in successful expansion and enrichment of CD56 positive, CD3 negative NK cells (>96%) with minimal T cell contamination.
- ✓ Using G-Rex6M plates, CD38 CAR-NK cells with CD38 KO expanded approximately 400 fold by D22.
- ✓ Extending cell expansion until D27 increased expansion 4-8 fold*
 ✓ 1,600-3,200 fold expansion*

by D27

Proprietary & Confidential 103



Process Continuity = **Product** Continuity



Figure 1. NK-cell expansion rate is comparable between Grex 6M, 100M, and 500M vessels.

Figure 2. NK-cell potency is comparable between Grex 6M, 100M, and 500M vessels.

ScaleReady,

G-Rex[®] for Donor-Derived Allo NK







Maximize fold expansion and minimize cell processing burden



Golubovskaya, V., Sienkiewicz, J., Sun, J., Zhang, S., Huang, Y., Zhou, H., Harto, H., Xu, S., Berahovich, R., & Wu, L. (2023). Car-NK cells generated with mrna-lnps kill tumor target cells in vitro and in vivo. *International Journal of Molecular Sciences*, 24(17), 13364. https://doi.org/10.3390/ijms241713364



Clinical Translation of Pluripotent Cell-Derived Off-the-Shelf Natural Killer Cell Cancer Immunotherapy

Frank Cichocki¹, Ryan Bjordahl², Svetlana Gaidarova², Paul Rogers², Raedun Clarke², Brian Groff², Stacey Moreno², Ramzey Abujarour², Megan Robinson², Greg Bonello², Tom Lee², Weijie Lan², Betsy Rezner², Stewart Abbot², Darin Sumstad¹, Bruce Blazar¹, Daniel Shoemaker², Scott Wolchko², Dan Kaufman³, David McKenna¹, Bahram Valamehr², Sarah Cooley¹ and Jeffrey S. Miller¹

¹University of Minnesota Cancer Center, Minneapolis, MN ² Fate Therapeutics, San Diego, CA ³University of California San Diego, CA







ScaleReady.

G-Rex[®] for Hematopoietic Stem Cell (HSC) therapy

ScaleReady,

Proprietary & Confidential 107

GREX Cityof Hope.

NOVEL, EFFICIENT SYSTEM FOR CULTURE OF TRANSDUCED HEMATOPOIETIC STEM PROGENITOR CELLS (HSPC) FOR VECTOR COPY NUMBER DETERMINATION

Monica Torres-Coronado¹, Ilse Arciga¹, Diana L. Browning¹, Xiu-Li Li¹, Agnes Gardner¹, Sergio Branciamore², John A. Zaia¹, Angelo A. Cardoso¹

Center for Gene Therapy ¹ and Department of Computational and Quantitative Medicine ² of City of Hope, Duarte, CA

Establishing Seeding Conditions to culture HSPC in G-Rex Plates (TDX with GFP LVV)



ScaleReady,



NOVEL, EFFICIENT SYSTEM FOR CULTURE OF TRANSDUCED HEMATOPOIETIC STEM PROGENITOR CELLS (HSPC) FOR VECTOR COPY NUMBER DETERMINATION

Monica Torres-Coronado¹, Ilse Arciga¹, Diana L. Browning¹, Xiu-Li Li¹, Agnes Gardner¹, Sergio Branciamore², John A. Zaia¹, Angelo A. Cardoso¹

Center for Gene Therapy ¹ and Department of Computational and Quantitative Medicine ² of City of Hope, Duarte, CA

Expansion Rates and VCN of TDX-HSPC Cumulative Expansion Rate (Fold Increase) 100000 120 UTDX Standard d1-d7 10000 0 Expansion Rate (Fold Increase) -O- TDX20 d7-d14 80 1000 UTDX G-Rex 100 0 TDX20 UTDX TDX20 UTDX d1 d7 d14 Standard G-Rex Days Post-TDX VCN by c-frag ddPCR UTDX ns Standard TDX20 ns Cell Viability % UTDX G-Rex 00 TDX20 000 <u>_</u> Standard G-Rex Standard G-Rex d14 d7 **Days Post-TDX** d14 Post-TDX d7 Post-TDX

- The Cell Viability of G-Rex cultures was consistently high during culture (d7, 96.8%; d14, 94.6%) and comparable to that seen in Standard cultures.
- Higher cell expansion rate was observed in G-Rex method during the first week, with similar expansion rates in the second week of culture.
- ✓ This novel, G-Rex-based system of TDX-HSPC expansion/differentiation for VCN determination is simple and robust, providing marked labor and cost savings.
- ✓ This system is being integrated in the QC testing for release of TDX-HSPC Drug Products for clinical manufacturing. This method is **suitable for automation** using a robotic platform integrating genomic DNA preparation and ddPCR processing.

ScaleReady.

G-Rex[®] for Pancreatic Islet Cells



Improved Islet Culture on Top of Gas-Permeable Membranes

Klearchos K Papas¹, Bernhard J Hering¹, Phillip R Rozak¹, John Wilson², and Efstathios S Avgoustiniatos¹ ¹Diabetes Institute for Immunology and Transplantation, Department of Surgery, University of Minnesota, Minneapolis, MN, USA ²Wilson Wolf Manufacturing, New Brighton, MN, USA



15 20 25 30

Time Posttransplant (Days)

10

ScaleReady.



Improved Islet Culture on Top of Gas-Permeable Membranes Klearchos K Papas¹, Bernhard J Hering¹, Phillip R Rozak¹, John Wilson², and Efstathios S Avgoustiniatos¹

¹Diabetes Institute for Immunology and Transplantation, Department of Surgery, University of Minnesota, Minneapolis, MN, USA ²Wilson Wolf Manufacturing, New Brighton, MN, USA



G-Rex[®] is an integrated development and manufacturing platform

ScaleReady.
G-REX: An Integrated Manufacturing Device



G-Rex[®] Integrated Development and Manufacturing Platform for CGT









✓ G-Rex is a "factory in a jar"

 Incorporate multiple unit operations within the G-Rex to streamline manufacturing process and reduce risk of contamination

Gagliardi C, Khalil M, Foster AE. Streamlined production of genetically modified T cells with activation, transduction and expansion in closed-system G-Rex bioreactors. Cytotherapy. 2019 Dec;21(12):1246-1257. DOI: 10.1016/j.jcyt.2019.10.006.

Joseph, Gil, et al. *Innovative Development of Closed CAR-T Platform.* WuXi Advanced Therapies. Accessed 12 Jan. 2021.

G-Rex[®] = Highly Controllable Process





CAR-T cell expansion platforms yield distinct T cell differentiation states

Hannah W. Song ● Michaela Prochazkova ● Lipei Shao ● ... David F. Stroncek ● Javed Khan ● Steven L. Highfill & 🖾 ● Show all authors

Open Access • Published: March 14, 2024 • DOI: https://doi.org/10.1016/j.jcyt.2024.03.003

- ✓ **G-Rex**[®] is the **most controllable** C> platform in the industry
- G-Rex[®] enables tight control over critical process parameters to deliver proper critical quality attributes

✓ G-Rex[®] is only platform that can be leveraged to deliver the proper cell quality in the proper cell quantity <u>at any scale</u>



Getting started in G-Rex®

ScaleReady,

Optimize expansion and read-outs at small-scale.



G-Rex "M" series well plates available in 2cm², 5cm², and 10cm²

ScaleReady,

G-Rex "M" Series Well Plate Experiment:

- Correlate glucose consumption and lactate production trends to cell density 1.
- Understand how culture duration affects the quality cell output 2.
- Confirm optimal minimum and maximum cell densities 3.
- Perfect cytokine supplementation strategy 4.

Linear Scale Up



G-Rex10M-CS

G-Rex50M-CS

G-Rex100M-CS



G-Rex500M-CS

122

G-Rex6M Well Plate® Optimization Study



Legend

Pink = Glucose [mmol/L] Blue = Lactate [mmol/L] Black = Cell Density (cells/cm²)

Circle Icon= 0.5x10⁶ cells/cm² seeding density

Triangle Icon = 1.0×10^6 cells/cm² seeding density

Square lcon = 2.0×10^{6} cells/cm² seeding density

ScaleReady,

 \equiv

G-Rex® vs. Static Gas Permeable Bags



Making Cancer History®



Figure 1. Assessment of TIL growth after rapid expansion in the traditional flask and bag vs. Gas-permeable flask (G-Rex) reveals a trend toward improved TIL expansion when using the G-Rex flask. (**A**) Fold expansion of post-REP TIL lines expanded shows a trend toward a better expansion when using the G-Rex. N = 10. (**B**) The G-Rex flask facilitates the expansion of TIL lines from which the growth is impaired in the REP using the traditional flask and bag devices N = 4. Statistical significance was determined by using a paired *t*-test.

					Bag-base	ed process	G-Rex-b	ased process
Item	Manufacturer	Catalog no.	Size	List price	Required per run	Cost per run	Required per run	Cost per run
TexMACS GMP medium	Miltenyi Biotec	170-075-306	2 L	\$ 305.00	1.00	\$ 305.00	1.00	\$ 305.00
MACS GMP Vectofusin-1	Miltenyi Biotec	170-076-165	1 mg	\$875.00	0.00	\$ -	1.00	\$ 875.00
G-Rex 100MCS	Wilson Wolf Corp	81100-CS	3 pack	\$729.84	0.00	\$ -	0.33	\$243.28
MACS GMP CD3 pure	Miltenyi Biotec	170-076-116	1 mg	\$ 1750.00	0.58	\$ 1015.00	0.08	\$ 140.00
MACS GMP CD28 pure	Miltenyi Biotec	170-076-117	0.5 mg	\$ 895.00	1.45	\$ 1297.75	0.20	\$ 179.00
hIL-7, premium grade	Miltenyi Biotec	130-095-363	$100 \ \mu g$	\$1150.00	0.13	\$ 149.50	0.45	\$ 517.50
hIL-15, premium grade	Miltenyi Biotec	130-095-765	$100 \ \mu g$	\$1150.00	0.04	\$ 46.00	0.15	\$ 172.50
Retronectin, GMP grade	Takara	T202	2.5 mL	\$1442.00	0.50	\$ 721.00	0.00	\$ -
Cell culture bags	Saint Gobain	290-C	10 pack	\$ 993.51	0.30	\$ 298.05	0.00	\$ -
Transduction bag	Saint Gobain	290-AC	10 pack	\$ 993.51	0.10	\$ 99.35	0.00	\$ -

The cost of generating 1.4×10^9 transduced cells in the bag- and G-Rex–based process. Cost for the bag process was calculated based on scaling up to a process that could generate as many cells as the G-Rex–based process.

Table 2. Bag and G-Rex process time comparison.

Operation	Bag-based process		G-Rex-based process		
	Method	Time (h)	Method	Time (h)	
Activation	Anti-CD3/CD28 in solution	1.00	Anti-CD3/CD28 in solution	1.00	
Transduction	Culture wash, transfer to Retronectin-coated bags	4.00	Vectofusin-1 in solution	0.50	
Transduction stop	Culture wash, transfer bag	2.00	Medium addition	0.50	
Expansion	Feed and transfer, as needed	4.00	NA	0.00	
Harvest	Centrifuge	2.00	GatheRex	0.50	
In-process monitoring	Mix cell suspension, sample for count (x6 d)	3.00	Sample supernatant for lactate (x6 d)	1.50	
		16.00		4.00	

The estimated time, in hours, required for each operation in the bag- and G-Rex-based process. Estimates are based on a single operator working in a non-GMP setting. NA, not applicable.

Gagliardi, C., Khalil, M., & Foster, A. E. (2019). Streamlined production of genetically modified T cells with activation, transduction and expansion in closed-system G-rex bioreactors. *Cytotherapy*, *21*(12), 1246–1257. https://doi.org/10.1016/j.jcyt.2019.10.006

ScaleReady,

Forget, M.-A., Haymaker, C., Dennison, J. B., Toth, C., Maiti, S., Fulbright, O. J., Cooper, L. J., Hwu, P., Radvanyi, L. G., & Bernatchez, C. (2015). The beneficial effects of a gas-permeable flask for expansion of tumor-infiltrating lymphocytes as reflected in their mitochondrial function and respiration capacity. *Oncolmmunology*, *5*(2). https://doi.org/10.1080/2162402x.2015.1057386

G-Rex[®] vs. Gas Permeable Bags





Figure 6. Large-scale G-Rex process compared with bag process lactate, transduction and expansion. (A, B) Lactate concentration (y-axis) plotted against process day (x-axis) for large-scale G-Rex100MCS (A) and bag cultures (B). (C–F) Comparison of transduction percentage (A), fold change (B), total cells at harvest (C) and total transduced cells at harvest (D) of cells generated in the G-Rex and bag-based processes. Each marker represents 1 donor. *** $P \le 0.001$, not significant P > 0.05.

Equivalent tdx %

- Significantly greater fold expansion in G-REX[®]
- Significantly greater total cell yield in G-REX[®]
- Significantly greater CAR+ cell yield in G-REX[®]

cant. Expansion in the G-Rex bioreactor resulted in significantly more total viable cells and transgenic cells at harvest. The average harvest from the G-Rex process was $1.4 \pm 0.1 \times 10^9$ transgenic cells. A hypothetical 80-kg patient who needs a dose of 5.0×10^6 cells/kg could receive >3 doses from one manufacturing process with the G-Rex. In comparison, a bag process, starting with the same number of cells, results in barely enough transgenic cells for a single dose of the same size. Further, the G-Rex-based process reduced the cost of materials by 38% and hand-on time by 75% to generate a batch of 1.4×10^9 transgenic cells.

ScaleReady.

G-Rex[®] for Donor-Derived Allo NK



Cytotherapy. 2012 October; 14(9): 1131-1143. doi:10.3109/14653249.2012.700767.

Large-scale *ex vivo* expansion and characterization of natural killer cells for clinical applications

NATALIA LAPTEVA¹, APRIL G. DURETT¹, JIALI SUN¹, LISA A. ROLLINS¹, LESLIE L. HUYE¹, JIAN FANG¹, VARADA DANDEKAR¹, ZHUYONG MEI¹, KIMBERLEY JACKSON¹, JUAN VERA¹, JUN ANDO¹, MINHTRAN C. NGO¹, ELAINE COUSTAN-SMITH³, DARIO CAMPANA³, SUSANN SZMANIA⁴, TARUN GARG⁴, AMBERLY MORENO-BOST⁴, FRITS VANRHEE⁴, ADRIAN P. GEE^{1,2}, and CLIONA M. ROONEY^{1,2}

¹Center for Cell and Gene Therapy, The Methodist Hospital, Texas Children's Hospital, Houston, TX

²Department of Pediatrics, Baylor College of Medicine, Houston, Texas, USA

³National University of Singapore, Singapore

⁴University of Arkansas Medical Center, Little Rock, Arkansas, USA

Figure 2.

G-Rex yields higher numbers of NK cells and requires fewer culture manipulations than

cultures in bags. (A) Schematic view of NK cell expansion in bags and G-Rexes. Less cell processing is required to grow a similar number of cells in G-Rexes. (B) NK cells from three different donors were seeded in the 197-mL bags and G-Rex100s. Bags were seeded with PBMC containing 2×10^6 NK cells and 2×10^7 K562-mb15-41BBL cells in a 40-mL starting volume and cultures required feedings with 40 mL fresh media and cytokine every other day, i.e. days 2, 4, 6 and 8, resulting in 200 mL of total volume by day 8. G-Rex100s were seeded with PBMC containing 2×10^6 NK cells and 2×10^7 K562-mb1L15-41BBL in 400 mL medium and were not subsequently fed. Cells were counted and analyzed by flow cytometry for percentage CD56⁺ CD3⁻ cells (n = 3). (C) Fold expansion of CD56⁺ CD3⁻ cells in G-Rexes and bags (n = 3). *P < 0.5.

G-Rex[®] vs. "Automated" Perfusion Bags

ScaleReady,

🜔 cytiva



! Unnecessary complexity





ScaleReady.



Day	Donor 1	Donor 2	Donor 3	Donor 4
5	6	6	6	6
6	12	6	6	6
7	12	8	6	6
8	12	8	8	6
9	12	8	8	6
10	12	8	10	6
11	12	8	10	6
12	12	8	12	6
13	12	8	12	6
14	12	8	12	6

Supplementary Table 4: Rocking rate settings (rocks per minute) used for Xuri W25 manufacturing process.

Supplementary Table 3: Media feed volumes and perfusion rate settings used for Xuri W25 manufacturing process. Innoc = inoculation.

Day	Donor 1	Donor 2	Donor 3	Donor 4
5	400 mL innoc.	400 mL innoc.	400 mL innoc.	400 mL innoc.
6	Add 600 mL, batch feed (=1 L total)	Add 600 mL, batch feed (=1 L total)	Add 600 mL, batch feed (=1 L total)	Add 350 mL, batch feed (=750 mL total)
7	1 L/day	0.5 L/day	1 L/day	Add 250 mL batch feed (=1L total)
8	1.5 L/day	0.5 L/day	1.5 L/day	no feed
9	1.5 L/day	0.75 L/day	1.5 L/day	0.25 L/day
10	2 L/day	0.75 L/day	2 L/day	0.25 L/day
11	2 L/day	1 L/day	2 L/day	0.25 L/day
12	2 L/day	1 L/day	2 L/day	0.25 L/day
13	2 L/day	1 L/day	2 L/day	0.25 L/day
14	Harvest	Harvest	Harvest	Harvest

Xuri cannot deliver consistency because of inherent design flaws

- Variability in rocking settings, feed volumes, and perfusion rates across donors
- Variability in rocking settings, feed volumes, and perfusion rates within donors

	Total ce	ll number (× 1	10 ⁹), Day 14	
Donor	Bag	G-Rex	Xuri	Prodigy
1	54.04	17.10	21.40	4.81
2	37.76	8.64	6.86	3.74
3	15.04	11.88	3.61	3.58
4	6.44	7.83	1.26	5.21
Avg	28.32	9.45	8.28	4.18
SD	21.65	4.21	9.04	0.90

ScaleReady,

G-Rex[®] vs. Automated Perfusion Bag



A SUPPLY CHAIN CRISIS STORY: CULTURE BAG SHORTAGE ENFORCED VALIDATION OF AN ALTERNATIVE EXPANSION SYSTEM FOR CAR T CELLS

Mercy Gohil^{1,2}, Don Hasenmayer^{1,2}, Kathleen Haines^{1,2}, Ankita Jain^{1,2}, Lauren Lewitt^{1,2}, Kevin Sporici^{1,2}, Edward C. Pequignot^{1,2}, Anlan Dai^{1,2}, Dmitri Negorev^{1,2}, Shane Mackey^{1,2}, Vanessa Gonzalez^{1,2}, Irina Kulikovskaya^{1,2}, Rachel Reynolds^{1,2}, Tyler Migliaccio^{1,2}, Andrea L Brennan^{1,2}, Theresa Colligon^{1,2}, Athena Russell^{1,2}, Ziming Wang^{1,2}, Carl H. June^{1,2,3}, Donald L. Siegel^{1,2,3}, Bruce L. Levine^{1,2,3}, Joseph A. Fraietta^{1,2,4}, Julie Jadlowsky^{1,2,3}, Gabriela Plesa^{1,2}, Megan M. Davis^{1,2,3}

Perelman School of Medicine at the University of Pennsylvania, Philadelphia, PA USA
 Center for Cellular Immunotherapies
 Department of Pathology and Laboratory Medicine
 Lepartment of Microbiology

Center for Cellular Immunotherapie mercy@upenn.edu

RESULTS

BACKGROUND

The recent supply chain crisis highlights a need to establish alternative manufacturing (MFG) protocols ensuring continuity of existing and new cell therapy (CT) clinical trials. Our mitigation strategy aimed at process redundancies to overcome creating production challenges due to limited availability of cell expansion culture bags for the Wave bioreactor, a critical unit of operation that we have used to successfully MFG thousands of gene-modified T cell products for 30+ clinical trials. Manufacturing Mesothelin (Meso) CAR T Cells from consented patient cryopreserved apheresis in the closed system 1L GREX required 1/4 of starting material, 1/5 of media and decreased manual effort through the culture duration compared to our standard WAVE manufacturing. We generated in the 1L GREX, strikingly 2-5 billion T cells from starting seed of 50 million elutriated lymphocytes, with sufficient transduction to meet most protocol-specified cell therapy doses and led to feasibility studies of replacing a full apheresis with peripheral blood collection as starting material.

GOALS

 Validate alternative ex vivo culture system for CAR T cells.

 Explore feasibility of using peripheral blood draw as starting material compared to full apheresis

METHODS

Meso CAR T cells were manufactured in parallel via the G-REX or conventional Wave bioreactor using consented patient starting material. In addition, Meso CAR T cells were expanded in parallel from healthy donor matched CD4+ CD8+ enriched cells from full apheresis or from peripheral blood draws. Critical quality attributes of the final T cell products, including viability, transduction efficiency, and phenotype were assessed. In vitro and in vivo functional studies are in progress.





 Image
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +
 +

 -</th

Figure 2. Phenotype of elutriated lymphocytes expanded in 1L GREX compared to WAVE.



Fig.2. Hierarchical charts depicting percentage of marker from flow cytometry OP 11 at A) 50 1 LGREX 80 10 00 MWE and 6) D0 Buseline. For P1 1 ools way and a start of the transformed and the transformed

Figure 3. CAR+ T Cell Expansion from Peripheral Blood as starting material compared to Apheresis



Fig.3. A) CD4+CD8+ T cells were enriched from healthy donor matched cryopreserved apheresis or PBMC. Enriched cells were seeded in LGREX, stimulated, transduced and expanded for 9 days. B) Table of mean± SD of fold expansion, % transduction and % CD3CD45 at 109 of apheresis and Peripheral blood starting material.

Conclusions

1L GREX required ¼ less seed number, 1/5 less media and 1/3 less consumables and had population doubling level 6.

CART cells with sufficient fold expansion and % transduction to meet protocol specified dose could be generated from GREX expansion of enriched T cells from healthy donor peripheral blood.

Acknowledgment We thank the Translational and Corelative Studies Laboratory, the Clinical Cell an Vaccine Production Facility, the Human Immunology Core, and all of our collaborator within the Center for Cellular Immunotherapies at the University of Pennstyvania, adm with the clinical full subjects who have conserved to the use of their samples for research

A SUPPLY CHAIN CRISIS STORY: CULTURE BAG SHORTAGE ENFORCED VALIDATION OF AN ALTERNATIVE EXPANSION SYSTEM FOR CAR T CELLS

M. Gohil¹, D. Hasenmayer¹, K. M. Haines¹, A. Jain¹, L. Lewitt¹,
K. Sporici¹, A. Dai¹, D. Negorev¹, V. E. Gonzalez¹, I. Kulikovskaya¹,
R. Reynolds¹, T. Migliaccio¹, A. L. Brennan¹, T. A. Colligon¹,
A. L. Russell¹, Z. Wang¹, C. June¹, D. Siegel¹, B. L. Levine¹, J. A. Fraietta¹,
J. Jadlowsky¹, G. Plesa¹, M. M. Davis¹
¹Center for Cellular Immunotherapies, University of Pennsylvania

Perelman School of Medicine, Philadelphia, PA, United States.

Using G-Rex[®] vs. WAVE...

- Reduces media consumption by 80%
- Reduces starting material requirements by 66%
- Reduces manual effort
- Reduces reliance on specialized capital equipment
- Reduces footprint
- Increases throughput (via enabling simultaneous manufacturing of products)

ScaleReady,

G-Rex[®] vs. Automated Perfusion Bag





G-Rex[®] vs. Automated Perfusion Bag



133



FIG. 4: Comparison of NK cell expansion in G-Rex flasks and WAVE bioreactor. PBMCs from healthy donor containing 5×10⁶ NK cells were co-cultured with K562-41BBL-mbIL-15 cells for 9 days. Cells were moved to WAVE on day 4 of culture. (A, B, C) Similar numbers and fold of NK expansion were observed in G-Rex and WAVE cultures. (D) Potency (killing K562 cells) of both products was also comparable. Similar expansion was obtained for three donors and the representative data for one donor is shown. Critical Reviews[™] in Oncogenesis, 19(1-2):121–132 (2014)

Clinical Grade Purification and Expansion of Natural Killer Cells

Natalia Lapteva^{1*}, Susann M. Szmania², Frits van Rhee², & Cliona M. Rooney^{1*}

¹ Center for Cell and Gene Therapy, Baylor College of Medicine, Houston, TX 77030; ²University of Arkansas Medical Center, Little Rock, AK 72205

toxicity assays). Interestingly, the WAVE-generated products contain fewer CD3⁺ T cells and a higher frequency of CD56⁺CD3⁻ NK cells, perhaps because T cells prefer static culture (Fig. 4). The clinical trial in multiple myeloma demands both autologous and allogeneic NK cells. In the case of allogeneic products, CD3⁺ T cells are efficiently removed (as well as CD3+56+ cells) by depletion with the CliniMACS CD3 reagent, reducing the CD3+ cell frequency from 19±9% to 0.4±0.5% (N=5). This depletion allows less than 5×10⁵ CD3⁺ CD56⁻ T cells kg⁻¹ for infusion, as mandated by the clinical protocol.

ScoleReady. Lapteva, N., Szmania, S. M., van Rhee, F., & Rooney, C. M. (2014). Clinical grade purification and expansion of Natural Killer Cells. Critical Reviews in Oncogenesis, 19(1–2), 121–132. https://doi.org/10.1615/critrevoncog.2014010931

Due to the fact that CD8+ effector cells are our desired cell population the increased representation of CD4+ cells in rapid expansions carried out in the WAVE bioreactor compared to static bags is problematic. The underlying cause of this skewed expansion could relate to differences in the culture microenvironment generated in the two styles of bioreactor. In static bags the cells rest on the lower surface of the bag and there is minimal mixing of the bag contents. In this situation there is extensive cell to cell contact and any secreted factors essential to cell expansion will exist as gradients, with the highest concentration surrounding the cells. Contrast this with the WAVE bioreactor, where the coupling of rocking with media perfusion creates a more homogenous suspension of cells and oxygenated nutrient rich media. In this situation any secreted factors will be diluted by the continual perfusion, while the constant motion will disrupt any gradient formation and reduce the overall time that cells spend in contact with each other. A second

There are some disadvantages with the WAVE bioreactor including a difficult transition from research scale expansions, where the initial cell number and seeding density are determined, along with the day of WAVE inoculation and rocking speed, to full scale clinical expansions. In addition to the WAVE bioreactors themselves it is necessary to purchase additional ancillary equipment. This, along with the inherent complexities of a system that depends on constant motion, such as a predisposition to electrical and mechanical failure, means that there is a need for multiple bioreactors to expand a single patient's cells. In order to exclusively adopt the WAVE bioreactor as the sole platform for rapid expansions, therefore requires a substantial initial investment, which may be beyond the means of a small cell production facility that produces only a minimal number of cell products; however, the WAVE bioreactor is ideal for larger scale manufacturing facilities, where multiple cell products need to manufactured concurrently. We have identified alternative bioreactors which don't require an additional expenditure on ancillary equipment, as they utilize equipment and monitoring systems that are standard to most cell production facilities [34].

G-Rex[®] vs. T-Flask

ScaleReady,

Proprietary & Confidential 135

 \equiv





FIGURE 6. Bioprocess optimization for CTLs production decreases the complexity, risk of contamination, and cost of manufacture. The conventional method of CTL culture using 24-well plates is convoluted, demanding frequent manipulation to sustain culture growth. This increases the manufacture cost and risk of contamination. In contrast, CTL culture using the optimized APC:CTL ratio in the G-Rex can produce the same number of cells in a shorter period of time thus decreasing the cost and complexity of manufacture.

The G-Rex Reduces the Time, Complexity, and Cost Associated With CTL Production for Clinical Use

The broader implementation of T-cell immunotherapy is limited by (i) the cost of production; (ii) the complexity of production, including repeat feeding of open culture systems, and multiple skilled "judgment calls," thereby limiting scalability, and (iii) the time required for activation and expansion. Figure 6 shows the effects of the G-Rex on these obstacles. The "hands on" time for CTL manufacture is reduced because fewer cell manipulations are required, which consequently reduces the labor costs, decreases the complexity of manufacture, and diminishes the risk of contamination. The duration of CTL manufacture is also reduced by the increased rate of expansion. Hence production of 1×10^{10} CTLs, which would typically require approximately 60 days and 129 hours of technician time by conventional culture, is reduced to 23 days with 3 hours of labor and a reduction in the number of interventions from > 200,000 to just 34, translating to a saving of >50% in CTL production costs (Fig. 6).

ScaleReady.

Vera, J. F., Brenner, L. J., Gerdemann, U., Ngo, M. C., Sili, U., Liu, H., Wilson, J., Dotti, G., Heslop, H. E., Leen, A. M., & Rooney, C. M. (2010). Accelerated production of antigenspecific T cells for preclinical and clinical applications using gas-permeable rapid expansion cultureware (G-Rex). Journal of Immunotherapy, 33(3), 305–315. https://doi.org/10.1097/cji.0b013e3181c0c3cb

G-Rex[®] vs. T-flask

\frown			- Sistem	a Socio Sanitario
GREX	386///	Ospedale di Bergamo	*	Regione Lombardia
			ASST	Papa Giovanni XXIII

Optimization of therapeutic T cell expansion in G-Rex device and applicability to large-scale production for clinical use

Elisa Gotti¹, Sarah Tettamanti², Silvia Zaninelli¹, Carolina Cuofano¹, Irene Cattaneo¹, Maria Caterina Rotiroti^{2,†}, Sabrina Cribioli³, Rachele Alzani⁴, Alessandro Rambaldi¹, Martino Introna^{1,*}, Josée Golay^{1,5}

¹ Center of Cellular Therapy "G. Lanzani," Division of Hematology, Azienda Socio Sanitaria Territoriale Papa Giovanni XXIII, Bergamo, Italy
 ² Tettamanti Research Center, Department of Pediatrics, University of Milano-Bicocca/Fondazione MBBM, Monza, Italy
 ³ Accelera srl, Nerviano, Italy

- ⁴ Nemine Medical Crimero Nemine
- ⁴ Nerviano Medical Sciences, Nerviano, Italy
- ⁵ Fondazione per la Ricerca Ospedale Maggiore, Bergamo, Italy



Figure 1. Expansion of CIKs from PB is more efficient in G-Rex devices compared with T-flasks. PBMCs were plated at 0.5×10^6 /cm² in G-Rex vessels (G-Rex-6M, -10M or -100M) or at 3×10^6 /mL in T-flasks in CIK conditions. After 10 to 11 days for G-Rex or 21 days for flasks, cells were collected and counted (A). The number of live CD3⁺ was measured by immunophenotyping and flow cytometry (B). Lactate levels in the supernatant of G-Rex cultures were measured at different times to follow cell growth. The results of six representative cultures are shown (C). (Color version of figure is available online.)

G-Rex[®] produces more cells in less time

We observed that <u>culture in G-Rex vessels in</u> these conditions <u>allowed more efficient</u> and rapid CIK <u>expansion compared with T-flasks</u>, starting from the same number of cells. Indeed, using 5×10^6 PBMCs, a mean of 320×10^6 total cells and 300×10^6 CD3⁺ T cells were reproducibly obtained in only 10 to 11 days of culture using the G-Rex 10M vessel (100 mL), compared with a mean of 140×10^6 total cells and 110×10^6 T cells obtained in 21 days in standard flasks (Fig. 1A and B). The rate of CIK expansion in G-Rex was quite reproducible, as shown by measured lactate levels over time in different batches (Fig. 1C) and cell yields on days 10 and 11 (Fig. 1A and B).

ScaleReady.

G-Rex[®] vs. "Automated" T-Flask







✓ Simplicity

1016 -

ScaleReady.

G-Rex[®] vs. Automated T-Flask



! Inefficient use of space! No high throughput manufacturing





✓ Efficient use of space
 ✓ High throughput manufacturing

ScaleReady,

Lonza





ScaleReady,

Lonza



➢ G-Rex[®] produces more cells than "automated" T-Flasks

G-Rex[®] produces cells with higher viability than "automated" T-Flasks

G-Rex[®] produces more CAR+ cells than "automated" T-Flasks

ScaleReady.

G-Rex[®] Throughput



Throughput In parallel processing with limited capital investment

Input Process Par	ameters
Manufacturing Lines (n)	1
Selection (hr)	2
Activation (hr)	24
Transduction (hr)	24
Expansion (d)	5

Drug Name	Est. Incident Rate (US)	Est. Eligible Patient Pop (US)
Carvykti®	35,000	14,000
Kymriah®	60,000	8,000
Yescarta®	60,000	8,000
Breyanzi®	60,000	8,000
Abecma®	35,000	5,000
Yescarta®	20,000	2,000
Breyanzi®	20,000	2,000
Tecartus®	5,000	1,000
Kymriah®	6,000	800
Tecartus®	1,500	300
Yescarta®	2,400	240

Miltenyi Biotec			Ori biotech	TERU	LONZ Cell & Gene	a	
		Process			Out	put	
Patients	Batch Length (hr)	Batch Length (d)	Total Time (hr)	Batch / Hr	Batch / Day	Batch / Week	Batch / Year
1	170	7.08	170	0.0059	0.141	0.988	51
2	170	7.08	340	0.0059	0.141	0.988	51
5	170	7.08	850	0.0059	0.141	0.988	51
10	170	7.08	1700	0.0059	0.141	0.988	51
100	170	7.08	17000	0.0059	0.141	0.988	51
1,000	170	7.08	170000	0.0059	0.141	0.988	51
10,000	170	7.08	1700000	0.0059	0.141	0.988	51
100,000	170	7.08	17000000	0.0059	0.141	0.988	51
GREX							
		Process	2		Out	put	
Patients	Batch Length (hr)	Process Batch Length (d)	Total Time (hr)	Batch / Hour	Out Batch / Day	put Batch / Week	Batch / Year
Patients 1	Batch Length (hr) 170	Process Batch Length (d) 7.08	Total Time (hr) 170	Batch / Hour 0.0059	Out Batch / Day 0.14	put Batch / Week 0.99	Batch / Year
Patients 1 2	Batch Length (hr) 170 170	Process Batch Length (d) 7.08 7.08	Total Time (hr) 170 172	Batch / Hour 0.0059 0.0116	Out Batch / Day 0.14 0.28	put Batch / Week 0.99 1.95	Batch / Year 51 102
Patients 1 2 5	Batch Length (hr) 170 170 170 170	Process Batch Length (d) 7.08 7.08 7.08	Total Time (hr) 170 172 178	Batch / Hour 0.0059 0.0116 0.028	Out Batch / Day 0.14 0.28 0.67	put Batch / Week 0.99 1.95 4.72	Batch / Year 51 102 245
Patients 1 2 5 10 10 10 10 10 10 10 10 10 10 10 10 10	Batch Length (hr) 170 170 170 170 170 170	Process Batch Length (d) 7.08 7.08 7.08 7.08	Total Time (hr) 170 172 178 188	Batch / Hour 0.0059 0.0116 0.028 0.053	Out Batch / Day 0.14 0.28 0.67 1.28	put Batch / Week 0.99 1.95 4.72 8.94	Batch / Year 51 102 245 465
Patients 1 2 5 10 10 100	Batch Length (hr) 170 170 170 170 170 170 170 170 170 170	Process Batch Length (d) 7.08 7.08 7.08 7.08 7.08	Total Time (hr) 170 172 178 188 368	Batch / Hour 0.0059 0.0116 0.028 0.053 0.272	Out Batch / Day 0.14 0.28 0.67 1.28 6.52	put Batch / Week 0.99 1.95 4.72 8.94 45.65	Batch / Year 51 102 245 465 2,374
Patients 1 2 5 10 10 100 1,000	Batch Length (hr) 170 170 170 170 170 170 170 170 170 170 170 170 170	Process Batch Length (d) 7.08 7.08 7.08 7.08 7.08 7.08 7.08	Total Time (hr) 170 172 178 188 368 2168	Batch / Hour 0.0059 0.0116 0.028 0.053 0.272 0.461	Out Batch / Day 0.14 0.28 0.67 1.28 6.52 11.07	put Batch / Week 0.99 1.95 4.72 8.94 45.65 77.49	Batch / Year 51 102 245 465 2,374 4,030
Patients 1 2 5 10 10 100 1,000 10,000	Batch Length (hr) 170 170 170 170 170 170 170 170 170 170 170 170 170 170 170 170 170	Process Batch Length (d) 7.08 7.08 7.08 7.08 7.08 7.08 7.08	Total Time (hr) 170 172 178 188 368 2168 20168	Batch / Hour 0.0059 0.0116 0.028 0.053 0.272 0.461 0.496	Out Batch / Day 0.14 0.28 0.67 1.28 6.52 11.07 11.90	put Batch / Week 0.99 1.95 4.72 8.94 45.65 77.49 83.30	Batch / Year 51 102 245 465 2,374 4,030 4,332

G-Rex[®] v. Prodigy v. Xuri

ScaleReady,



✓ G-Rex[®] produces greater fold expansion in 4/4 donors when compared with Prodigy/Xuri







✓ G-Rex[®] achieves consistently higher cell concentrations than Prodigy or Xuri

- ✓ G-Rex[®] has greater efficiency with media/cytokines consumption
- ✓ **G-Rex**[®] lowers the cost of goods sold by reducing inputs and maximizing outputs



*Proper use of G-REX will yield 4x10⁶ cells/mL regardless of donor/patient. The authors intentionally used fed-batch or Prodigy-like processes in G-REX to maintain "consistency" in media feeding/exchange frequencies across platforms. Even with a suboptimal use of G-REX, the superiority over Xuri in terms of the efficiency of media usage is evident.



G-Rex[®] vs. Prodigy vs. Xuri





*G-Rex space utilization assumes 2 VIOS250i incubators loaded with 10 G-Rex500M-CS (each) for a total of 20 devices. Footprint is calculated based on surface area of incubators as reported by Thermo Fisher.

 ${}^{*}\!\mathsf{Xuri}\ \mathsf{space}\ \mathsf{utilization}\ \mathsf{is}\ \mathsf{calculated}\ \mathsf{based}\ \mathsf{on}\ \mathsf{rocker/pump}\ \mathsf{width}$

*Prodigy space utilization is cacluated based on Prodigy unit only (not extra space for tube welder/sealer & barcode scanner)









ScaleReady,

G-Rex[®] vs. Prodigy vs. Xuri





 ** High throughput capacity for GREX was calculated based on filling a VIOS250i incubator with maximum capacity of 45 G-Rex100M-CS (each, 90 total in stacked incubators)

**Prodigy/Xuri calculations assume 1 patient dose per instrument



ScaleReady.

Hannah W. Song, Michaela Prochazkova, Lipei Shao, Roshini Traynor, Sarah Underwood, Mary Black, Vicki Fellowes, Rongye Shi, Marie Pouzolles, Hsien-Chao Chou, Adam T. Cheuk, Naomi Taylor, Ping Jin, Robert P. Somerville, David F. Stroncek, Javed Khan, Steven L. Highfill, CAR-T cell expansion platforms yield distinct T cell differentiation states, *Cytotherapy* (2024), doi; https://doi.org/10.1016/j.jcyt.2024.03.003

Proprietary & Confidential 151

G-Rex[®] vs. Prodigy vs. Xuri





Supplementary Figure S5: Naïve/stem central memory-like and effector memory phenotype over time



 ✓ G-Rex[®] produces CAR-T product with equivalent % T_{n/scm} in normoxia conditions

 ✓ G-Rex[®] produces CAR-T product with greater % T_{n/scm} in hypoxia conditions and greater fold expansion

Figure 5: Hypoxia culture yields a more naïve/stem central memory ($T_{n/scm}$)-like and less activated phenotype, with longer time in culture corresponding to greater proportion of $T_{n/scm}$ -like cells. (A) Dissolved oxygen (DO) concentration in Prodigy,

ScaleReady,

G-Rex® Grant Program for the advancement of CGT development and manufacturing
CGT Investment Overview



CGT Investment Overview

nature biotechnology

Biotech financing: darkest before the dawn



ScaleReady,

Proprietary & Confidential



Cell Therapy Becaccon byhansonwade

Deal and Company Landscape



Deals by Year & Potential Total Deal Value

What do investors want to see?



What do investors want to see?



Figure 11. (A) Rank of factors selected as first-place for most influential in decision-making to invest in a CGT product development opportunity. (B) Rank of factors selected in the top three for most influential in decision-making to invest in a CGT product development opportunity. (Color version of figure is available online.)

Kunze-Küllmer, M., Goonewardene, A., Kili, S., Theoharis, S., & Rivers, P. (2024). Cell and Gene Therapy Investment: Evolution and Future Outlook on Investor Perspectives. Cytotherapy. https://doi.org/10.1016/j.jcyt.2024.02.017

What do investors want to see?



Figure 12. Rank of factors selected in the top three for greatest barriers in decision-making to invest in a CGT product development opportunity.

Kunze-Küllmer, M., Goonewardene, A., Kili, S., Theoharis, S., & Rivers, P. (2024). Cell and Gene Therapy Investment: Evolution and Future Outlook on Investor Perspectives. Cytotherapy. https://doi.org/10.1016/j.jcyt.2024.02.017

ScaleReady.

Proprietary & Confidential 160

Company A

- Excessive cash burn rates
 - Imprudent decision making
 - Inefficient use of people, space, and materials
- **!** Slow to generate clinically significant data
- ! Doesn't keep the main thing the main thing

Company B

- ✓ Modest and reasonable cash burn rate
 - ✓ Prudent decision making
 - Highly efficient use of people, space, and materials
- ✓ Quickly generates clinical significant data
- ✓ Keeps the main thing the main thing!



Where do we go from here?



Where do we go from here?

CGT needs to transition away from "all-in-one" bioreactors













ScaleReady,

3

Development

- X No native small-scale model
- ! Unpredictable translation of small-scale data
- Excessive time and money spent attempting to create poorly predictive scale down models
- Excessive time and money spent optimizing CMC at full-scale

WHY REINVENT THE WHEEL WHEN YOU DON'T HAVE TO?



Manufacturing

- X Continuous fluidic connection of drug product
- ! Constraint is the manufacturing process itself
- Excessive investment in capital equipment for minimal throughput
- Limited potential for further innovation



CGT needs to transition away from "all-in-one" bioreactors













CGT needs to embrace **simplicity** and recognize it as the ultimate **sophistication**



Development

- ✓ Native small-scale model
- ! Predictable translation
- ! Linear path of least resistance to IND
- ! Cost- and time-efficient CMC improvements

Manufacturing

- ✓ De-couple drug product from instruments
- ! Constraints are individual unit operations
- ! Modest investment for high throughput
- ! Future-proofed



Surface Area





G-Rex® Grant Program

for the advancement of CGT development and manufacturing





G-Rex[®] Grant Program: Who? / What? / Where? / Why?





Estimated Time Requirements *over 18 months* = 6 hours

G-Rex[®] Grant Partners





G-Rex[®] is the most future-proofed manufacturing platform in the industry





G-Rex[®] with integrated robotic automation

Author Email: asmith@markertherapeutics.com

Robotic Automation of T Cell Generation for the Treatment of Acute Myeloid Leukemia (AML)



Anastasiya Smith², John d'Aigle¹, Tara Shahim², Sherket Peterson¹, Thorsten Demberg², Derian Salas¹, Sam Chang², Aditya Tandon¹, Jeannette Crisostomo², Alex Riley¹, Tsvetelina Hoang², Jaime Avalos¹, Juan Vera², Jose-Manuel Collados¹. ¹Division of Healthcare, Consumer Segments & Service Robotics, ABB Inc, Houston, Texas 77021, USA; ²Marker Therapeutics, Inc., Houston, Texas, 77027, USA.



ABB Robotics for Healthcare

Robotic Automation of T-Cell Generation for the treatment of Acute Myeloid Leukemia (AML)

embedded video; does not work as PDF

G-Rex[®] & Robotic Automation



Cellular Origins Partners With ScaleReady to Simplify, Standardise, and Automate Cell Therapy Manufacturing

- Partnership aims to integrate ScaleReady's CGT manufacturing workflow with Cellular Origins' robotic technology
- Initially focussed on automation of ScaleReady's G-Rex platform

May 11, 2023 07:28 AM Eastern Daylight Time











- ✓ G-Rex[®] is being automated in ways common to other mature industries
- ✓ Robots can be trained to maneuver and manipulate G-Rex[®] with consistency
- ✓ "Lights off" manufacturing with a G-Rex[®] centric approach
- ✓ G-Rex[®] based approaches are agonstic to centralized or de-centralized manufacturing





- ✓ G-Rex[®] is being automated in ways common to other mature industries
- ✓ Robots can be trained to maneuver and manipulate G-Rex[®] with consistency
- ✓ "Lights off" manufacturing with a G-Rex[®] centric approach
- ✓ G-Rex[®] based approaches are agonstic to centralized or de-centralized manufacturing







Here, we describe the performance of a new prototype robotic system designed to operate industrystandard cell therapy manufacturing equipment. We robotically automate the culture of CD8+ T cells using standard G-Rex 10M CS and Xuri Cell Expansion System W25 platforms, in combination with a modified Heracell Vios 250i cell culture incubator. We show that cell yields, phenotype, and viability are comparable between robotic and manual processes. The use of industry-standard equipment in a robotic process can greatly increase the pace of both process development and scale up for novel cell therapy products.



Melocchi, A., Schmittlein, B., Jones, A. L., Ainane, Y., Rizvi, A., Chan, D., Dickey, E., Pool, K., Harsono, K., Szymkiewicz, D., Scarfogliero, U., Bhatia, V., Sivanantham, A., Kreciglowa, N., Hunter, A., Gomez, M., Tanner, A., Uboldi, M., Batish, A., ... Esensten, J. H. (2023). Development of a Robotic Cluster for Automated and Scalable Cell Therapy Manufacturing. https://doi.org/10.1101/2023.12.21.572854

